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
Fig. 2

### I - 1S142 KIT INCLUDING :

- 1 guidance cone
- 1 linking device
- 8 linking screws

Fig. 3


Fig. 4

### II - ASSEMBLING :

1) Smooth distal part of the positive (fig. 1) and fasten the 1S112 tip with 2 little nails (fig. 2), then shape again its initial contour using plaster. Smooth the assembly.

2) Screw on the direction rod in 1S112 (fig. 3) and place the 1S129 tool (fig. 4).

**CAUTION !** Take care of the place of the metal centring pins A. They will maintain the 1K90 cup and determine the position of the releasing fork (draw a mark).

3) Slide the silicone band B on the 1S129 tool C after isolation with release agent (fig.5). Take care of the mounting direction for this silicone band.

4) On 1S129 tool place the 1K90 cup that is centered by the 2 metal centring pins (fig. 6). Then stick the self-adhesive protecting chip on the bottom of this cup.

5) Laminate or thermoform the socket as usually (fig. 7).

6) Remove the socket from the plaster and dismantle the forming tool and the tip. In the socket, at 38 mm from the lower extremity of the 1K90 cup, bore a Ø 24 mm hole in front of the fork (fig.8) in order to be able to remove it. Remove the fork from 1S140 (drawing it with pliers).

7) Screw the quick fastening device 1S140 onto the linking device with the 4 supplied screws (fig. 9).


8) Mount the assembly inside the socket and fasten it to the 1K90 cup (supplied screws).

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