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
6C99
INSTRUCTIONS FOR USE

A CONTENT:


- 1) A set of two complete joints with uprights that can be positioned on the right side and on the left side and are fitted with a 0° extension stop (blue) and a 20° flexion stop (blue).
- 2) A small bag including :

<ul style="list-style-type: none"> 2 x 6° extension stops (beige) 2 x 10° extension stops (green) 2 x 15° extension stops (black) 2 x 20° extension stops (yellow) 2 x 30° extension stops (red) 2 x 40° extension stops (brown) 	and	<ul style="list-style-type: none"> 2 x 40° flexion stops (brown) 2 x 60° flexion stops (red) 2 x 80° flexion stops (yellow) 2 x 90° flexion stops (black) 2 x 100° flexion stops (green)
--	-----	---


B REPLACEMENT OF AN EXTENSION STOP

- 1) Unscrew and remove screw A
- 2) Remove the extension stop that must be changed
- 3) Replace it with the convenient extension stop
- 4) Re-screw screw A.

C REPLACEMENT OF A FLEXION STOP

- 1) Unscrew and remove screw B
- 2) Remove the flexion stop that must be changed
- 3) Replace it with the convenient flexion stop
- 4) Re-screw screw B.

D SOFTENING THE SPLINTS FOR FORMING

- 1) Tip : To ensure a correct alignment of the splints when making the orthosis, we advise you to use our 6C96 forming kit.
- 2) Softening in the oven : Set the oven to 140°C and heat the splints for 20 minutes. As the metal screw might damage the plastic, it is essential to remove each of them before heating the splints.
- 3) Softening with hot air blower : Do not heat to more than 150°C.
WARNING : Stop heating immediately if you see that the colour of the material is changing.

E COOLING

In all cases leave the splints to cool on the plaster cast without accelerating the cooling process with compressed air.

F CUTTING TO LENGTH - DRILLING

All the operations to prepare the splints (cutting to length, drilling, etc...) should be carried out only when the splints have completely cooled.

E ACCESSORIES

- 1) PROTEOR 6C96 forming kit : for easy trimming when forming.
- 2) 6C97 lock : This accessory allows for an easy and effective locking in any predetermined position.
- 3) 6C98 condyle cushions : Comfort is improved by this set of two condyle cushions with accessories for fastening to the splints

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