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- 1 Clapet anti-retour (x1)
- 2 Raccord coudé (x2)
- 3 Raccord droit (x2)
- 4 Tube silicone Ø2 x Ø6 (x1)
- 5 Tube polyamide Ø4 x Ø6 (x1)

## B) DESTINATION

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## B) USE

1S160 anti-return valve should be used only to manufacture lower limb prostheses.

## C) DESCRIPTION / INDICATION

Used with a liner and a sealing device –such as 1S500-, the anti-return valve pushes the air out of the socket during stance phase and prevents it from returning into the socket during swing phase. The resulting vacuum reduces shearing and friction forces, and therefore allows a better contact between the liner and the socket.

## D) IMPLEMENTATION

Laminate the socket with your standard procedure. At the location chosen to connect the anti-return valve, drill and tap a 4 mm Ø M5 hole. Screw the selected connector in this hole and check that it does not protrude into the socket. If necessary the connector length can be shortened with a cutter (fig. 1). Apply adhesive resin (PC232-5) on the thread of the connector, then screw it in the tapped hole of the socket. Connect silicone tube **4** to the connector and cut it at the desired length. Connect the anti-return valve **1** to the silicone tube, taking care of the air transfer direction (arrow, fig. 2).

**Particular case of a dick socket :** If the socket is very dick (such as a double socket), it is better to proceed as follows to ensure tightness :

- Drill a 6.2 mm Ø hole at the selected location
- Insert polyamide tube **5** in the hole and cut it to the right length
- Remove the tube, then apply adhesive resin (PC232-5) to it
- Re-insert the tube in the hole, let resin dry, then tap (M5) the inner diameter of the tube (from the outside face of the socket)
- Apply PC232-5 adhesive resin to the thread of the selected connector
- Screw it in the tapped hole previously drilled in tube **5**
- Finish assembling as indicated previously.

**E) NOTE** Delivered with two sets of threaded connectors, the anti-return valve can be used with both trial socket and final socket.

**F) MAINTENANCE** To operate well, the anti-return valve should be protected against dust. Regularly check that it operates well and replace it if necessary.

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