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



INSTRUCTIONS FOR USE 1S101

1S10199-12/08

I.) PREPARATION FOR THERMOFORMING OR LAMINATION (See Fig. 1)

- Fix Delrin template **1** on plaster cast **2** with the four supplied nails **3**.
- With a drilling machine fitted with a Ø2 drill (2.4mm) drill a hole using the central hole of Delrin template **1** as guide. To provide optimal vacuum, this hole must extend to the base of plaster cast **2**.
- Fix bottom **4** on Delrin template **1** with cross-head screw **6**. Check that tip **5**, which is inserted in the outlet hole **4.1** of bottom **4**, is so orientated that it allows the string to run with minimal friction against the socket relief.

II.) MAKING THE SOCKET

Thermoform or laminate as usual, taking care to carefully apply socket material around bottom **4**.

III.) SOCKET FINISHING (See Fig. 1 and 2)


After socket thermoforming or lamination :

- Cut off material of socket **7** around screw **6** and remove it.
- Break plaster and remove Delrin template **1** from the inside of the socket.
- Also cut off material of socket **7** at the end of tip **5**, then remove tip **5** to ensure free access to outlet hole **4.1** of bottom **4**.

IV.) STRING INSTALLATION (See Fig. 2 and 3)

- Thread string **8** through bottom **4** and let it go out through outlet hole **4.1**.
 - Find out the best position for string stopping device **9**
 - . positions **a** and **a'** are wrong
 - . position **b** is correct.
 - Fix string stopping device **9** on socket **7** with the two Ø3mm rivets supplied (use string stopping device **9** as a drilling template for the two Ø3mm holes drilled for the rivets).
 - Cut sheath **10** of string **8** to the desired length.
 - Thread string **8** in sheath **10**.
 - Insert sheath **10** in outlet hole **4.1** of bottom **4**.
 - Stick sheath **10** to the socket with Loctite 425.
- Important:** Knot the end of string **8** to prevent it from falling into sheath **10** by mistake.

V.) OPERATION (See Fig. 3 and 4)

Position  = free string Position  = locked string

Note: To make traction easier thanks to a better handling of string **8** when the socket is worn, handle **11** (included in the package) can be installed at its end.

VI.) USE / ADVICE

This product is designed for patients with high activity level. We recommend you to check regularly string wear.

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