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1S145

Rack Lock for Liners and Lamination Anchor

Technical Instructions

1 - Use

The 1S145 lock is designed to secure an appropriate liner in the socket. This lock is laminated in the distal part of the socket.

This product is designed exclusively for lower limb prosthesis.

The use of this lock is limited to patients weighing 100 kg or less. Carried loads and patient activity must be taken into account when calculating the weight.

2 - Description

The 1S145 lock comes with :

- 2 positioning dummies
- 1 lamination anchor M36x1.5
- 2 locking pins M10 (1 long and 1 short)
- 1 bending spacer

3 - Technical data

- Products recommended for assembly

Lower connection : Adapter of 1K182 or 1K183 type
Threaded union M36.

- Authorized environmental conditions

Use temperature : -20°C à +65°C.

Relative air humidity : 0% to 90%, no condensation.

- Prohibited environmental conditions

Mechanical vibrations or shocks.

Sweating, urine, fresh water, seawater, acids, alkalis.

Dust, sand, strongly hygroscopic particles (eg : talc).

3-1 - Product life

This component has been tested in accordance with standard ISO10328, load level P5 (100 Kg / 220 lbs) ; the cyclical test of 3 million cycles corresponds to a use of 4 or 5 years depending on the patient activity.

3-2 - Warnings

In case of assembly with unauthorized prosthetic components:

- Risk of injury due to product breakage or warping

- Use only prosthetic components that have been approved for use with this product.

Read the instructions for use of the prosthetic components to check that they can be assembled with each other.

In case of use in unauthorized environmental conditions:

- Risk of injury due to product damage.

- Risk of malfunction

- If the product has been exposed to unauthorized environmental conditions, check that it is not damaged.

- In case of obvious or suspected damage, stop using the product.

- If necessary take appropriate measures (eg : cleaning, repair, replacement, control by the manufacturer or by a specialized workshop, etc ...).

4 - Final Assembly

- Position of the lamination anchor (3) : lateral or medial. If positioning is wrong, the lamination anchor (3) will have no support. This may cause a very high stress on the screwed adapter and break it. If the body of the locking mechanism (4) is too loose from the lamination anchor (3), it may break with in a risk of falling for the patient and thus injury.

- During alignment the body of the locking mechanism (4) must not be unscrewed from the lamination anchor (3) by more than half a turn.

- For a correct position of the lamination anchor (3), refer to figure 8. Make sure that the release button (7) is positioned on the side of the clamping screw of the lamination anchor (3).

4.1 - Assembly of the lamination anchor (3)

Shape the prongs of the lamination anchor (3) according to the socket outline.

- Bend the prongs of the lamination anchor (3) using the appropriate craft processes. In the event of improper implementation, the steel may crack.

- Never bend the prongs of the lamination anchor (3) when the body of the aluminum locking mechanism (4) is screwed because this may damage the aluminum thread and destroy it. When shaping the prongs of the lamination anchor (3), use either a steel or titanium adapter, or the 1K17910 steel bending spacer (6) screwed in M36 and tighten the cylinder head screw at 10 Nm. After bending, the area of the prongs that has been shaped must be polished to remove rough edges.

-The lock must always be placed in the extension of the stump, never in the buid line. This may require assembling a translation adapter.

Remove the locking mechanism (4) from the lamination anchor (3). Coat the thread of the lamination anchor (3) and of the cylinder head screw with a silicone release agent, then screw the positioning tool (2).

During lamination, the positioning tool (2) must be screwed and the cylinder head screw must be tightened to 4 Nm. Protect the slot of the lamination anchor with Plastiline.

If you have placed a rod in the cast to position the lock, place the spherical part of the positioning tool (1) on this rod in the plaster positive. If you are working without a positioning rod, attach the spherical part of the positioning tool (1) to the appropriate position. Insulate the plaster and spherical part of the positioning tool (1) with a PVA bag. All reinforcement materials placed before the lamination anchor (3) must be bound to the base of the spherical positioning tool (1). Place the lamination anchor (3) with the positioning tool (2) screwed onto the spherical positioning tool (1).

Make sure that the direction of the lamination anchor (3) has not changed.

The reinforcement materials which have been placed after the lamination anchor (3) can be bound directly to the positioning tool (2). Insulate the reinforcement with a PVA bag and laminate the socket.

Remove positioning tool (2), plaster and spherical positioning tool (1). Also remove the laminated resin edge at the distal part of the lamination anchor (3). Screw the body of the locking mechanism (4) until it stops in the lamination anchor (3) and tighten the clamping screw of the lamination anchor (3).

Tighten the cylindrical head screw with an appropriate wrench.

Tightening torque : 10 Nm

When finishing the prosthesis, coat the threads of the screw with Loctite® XC047. Place the adapter in the distal screw connection of the locking mechanism body (4). Check the function of the locking mechanism (4) with the appropriate locking pin (5).

4.2 - Definitive Assembly

Screw the 1K182 or 1K183 pyramid adapter fully into the M36x1.5 thread of the body.
Adjust orientation by unscrewing no more than a quarter of a turn.

Caution : Never unscrew the anchor more than a quarter of a turn and check that at least 7 threads are engaged.

Otherwise, the component may be weakened.

The distal screw connection between the lock and the adapter can be done with the M5 screw for fittings at the orthopedic workshop only. On delivery to the customer, the distal screw must always be secured with Loctite® XC 047.

Screw the locking mechanism with an appropriate wrench.

Tightening torque : 10 Nm

4.3 - Use

The liner and the lock are connected by the locking pin (5) that engages in the lock and holds the liner. Advise your patient that a correct alignment of the locking pin (5) with the locking mechanism (4) is very important. The locking pin (5) must not be forced into the locking mechanism (4). If the locking pin (5) is not correctly aligned with the locking mechanism (4), it may get stuck in the locking mechanism (4) and it will be difficult to release it. The operation of the locking mechanism (4) must be checked at regular intervals that do not exceed one year.

Placing the Liner

1. Check that the locking pin (5) and the liner are not damaged.
2. Tips: Align the locking pin (5) with the longitudinal axis of the stump to prevent damage. Roll the liner and place it on the distal end of the stump.
3. Unroll the liner on the stump, with no crease, no displacement of soft tissues and no air pocket.
4. Check the adjustment and direction of the liner.

Placing the socket

1. Put the socket on the liner until the locking pin (5) slides into the opening of the lock.
2. Insert the locking pin (5) completely into the lock.
3. Before using the prosthesis, check that the locking pin (5) has reached its stop position in the lock.

Taking off the socket

Press the release button (7) of the notched device and remove the stump and the liner from the socket.

5 - Maintenance

CAUTION : if the maintenance instructions are not respected, there is a risk of injury due to the modification or the loss of functions and to the damage of the product.

Please respect the following maintenance advice :

- After the first 30 days of use, visually inspect the components of the prosthesis and check their functions.
- The wear of the complete prosthesis should be checked during normal consultations.
- A safety check should be performed each year.

CAUTION : Check that nothing comes into the mechanism (dust, threads of a cup ...)

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