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PHONE : +33 3 80 78 42 10 - FAX : +33 3 80 78 42 15

www.proteor.com - htorthopedie@proteor.com

FOUR-AXIS KNEE JOINT WITH LOCK 1M102V Instructions for use

A) USE

This knee joint is designed for patients :

- who are moderately active and whose weight added to any load being carried does not exceed 100 kg
- who can progressively learn to walk " freely " during rehabilitation.
- who need an improved stability due to their activities (walking on very rough grounds,...).

Caution : This is not a lock knee. Only use the lock in the conditions indicated above.

The lock can be removed by a CPO if the patient does not use its function any more (Please see § K : TRANSFORMATION IN A FREE KNEE).

B) CONNECTIONS (Fig.1)

This knee can be assembled proximally (1) to any link piece for PROTEOR standard knee (1K03, 1K30, 1K160,...) and distally (2) to any Ø30 mm tube for prosthesis (1G01,...).


FIG 1

C) ALIGNMENT (Fig.1)

We recommend the alignment shown opposite.

D) EXTENSION ASSIST ADJUSTMENT (Fig.2)

Tighten or loosen screw 5 with a flat blade screwdriver to achieve the desired extension force.

E) EXTENSION BEARING ADJUSTMENT (Fig.2)

Caution : This adjustment affects the knee geometry and therefore the alignment. Undoing screw 6 improves flexion safety, but makes knee flexion harder. Make that adjustment with a 3 mm hex wrench when the patient weight is applied to the prosthesis in order to keep the bearing correctly in its housing.

F) JOINT FRICTION ADJUSTMENT (Fig.2 and 2-1)

Adjust only the two screws of the **right side** of the knee. Any adjustment made on the left side of the knee could damage the knee function.

- a) Unlock screw 3 with a 2 mm hex wrench.
- b) Screw in or out the main screw 4 with a 4 mm hex wrench to achieve the desired adjustment
- c) Then screw in screw 3 with a 2 mm hex wrench (tightening torque 2,5 N.m)


FIG 2


FIG 2.1

G) RECOMMENDATIONS

To eliminate noise in the cosmetic foam cover, do not use talcum, because it could damage the knee function and reduce the product life. Use preferably a silicone spray.

H) MAINTENANCE

Only the extension assist (ref : 1M10270) can be replaced.

CAUTION

In maximum flexion the socket must necessarily get in contact with the tube-holder.

Any other contact area is forbidden.
Check that the lock is not used as a bearing

WRONG




CORRECT

I) HANDLE

Use PROTEOR handle (ref : 1X110) or any other handle with at least a 12 mm stroke.

J) LOCK ADJUSTMENT

- a) Position the sheath and the cable, so that the latter be centered and aligned relative to the lock (Fig.3)
- b) Lock the knee and make sure that the lock is in contact with the linkage. In this position, the cable should be slightly slack (Fig.4). Using a 2 mm hex wrench, screw in the screw to achieve a 0,8 Nm torque.
- c) Adjust the lock position to have a 1,5 mm clearance between the lock and the linkage (See Fig. 5) by setting screw A with a 2,5 mm hex wrench (See Fig.3).

Before these adjustments, check that the knee is in extension.

Check that the knee is always conform to Fig. 5. Otherwise, an increased slack may damage the lock and a too large slack may make the patient fall.




In unlocked position, check that the lock is fully disengaged to allow flexion.

If the prosthesis alignment is changed, modify the lock engagement.

The knee should be locked in full extension. Otherwise the lock friction may damage the rear linkage. This device is not designed to get unlocked when the patient weight is applied to the prosthesis.

K) TRANSFORMATION IN A FREE KNEE

Caution : This operation is irreversible

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