# GENOU MONOAXIAL À FREIN 1M115 / 1M149

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
## Instructions for Use

### 1 – Description and functions

This knee is intended exclusively for lower limb prosthesis.

It is a single axis knee with hydraulic braking in stance phase, which provides 3 different adjustments and extension assist.

During stance phase, the ground reaction force closes a jaw that locks the rotation axis.

The following 3 adjustments must be performed :

**R** - Adjustment of the extension assist force

**F** - Adjustment of the joint friction value

**S** - Adjustment of brake sensivity according to load.

### 2 – Indications and maximum weight

It is specifically recommended for patients with medium to low activity level (L1/L2) up to 100 Kg (220 lbs) (carried load included).

This knee weighs 540 g for 1M115 model and 368 g for 1M149 model (titanium).

Maximum knee flexion : 150°. However it can be limited by the socket dimensions or by the cosmetic cover.

**CAUTION !**

In maximal flexion the socket must necessarily get in contact with the collar of tube-holder 3. Any other contact area is forbidden (see Fig. 1).

### 3 – Parts for connection and finishing (See Fig. 1)

RELATED PRODUCTS RECOMMENDED :

Upper connection with female pyramid adapter of 1K172 type.

Lower connection with adapter of 1D41 type.

FINISHING : One-part cosmetic cover 1G13 (complete or partial)

It is recommended to reinforce the cosmetic cover in front of the joint to increase its product life.

The cosmetic cover can modify the initial adjustments of the knee and limit some motions. In any case you should place the cosmetic cover before checking the knee functions.

### 4 – Alignments (See Fig. 2)



Function of braking mechanism essentially depends on prosthesis alignment. In the sagittal plane, the alignment shown opposite is highly recommended.

In the frontal plane the load line will come through the knee centre and the foot centre.

During normal walk after installation of the extension assist and after the different adjustments described below, the knee must be in extension during stance phase.

**Caution !** The load line must go in front of the rear axis (fig. 2) ; otherwise it could damage the knee prematurely and lead to a risk of breakage and therefore of falling.

## 5 - Installation of the extension assist (fig.3)

Install the parts of the extension assist as shown opposite : cable (a), adjustable bearing (R), mobile support (c) and springs (d).

**Caution !** If there is rubbing between cable (a) and lateral screws (2.1) (2.2) of adapter (2) during knee flexion, replace these screws by the shorter ones (12 mm long) supplied with the knee.

## 6 - Adjustments

This knee allows 3 possible adjustments.

To ensure the patient's safety, the prosthetist should perform the first fitting trials between parallel bars and with factory settings.

Then make the adjustments preferably in the following order:

1. Adjustment R of the extension assist force (fig. 3)
2. Adjustment F of the joint friction value (fig. 4)
3. Adjustment S of braking sensitivity according to load (fig. 5).

The foot model used can have an impact on the knee adjustments. The knee should be adjusted each time the foot is changed. The use of energy return feet such as Dynastar or DynaCity is recommended.

Ensure the patient's safety in all circumstances and especially during slope or stair descent.

A new adjustment of the knee is necessary after each maintenance operation.

### 1. Adjustment of the extension assist (fig.3)

Adjust the extension assist force by screwing or unscrewing adjustable bearing (R) with a 4 mm male hex wrench (k) :

- clockwise : tension will increase
- anti-clockwise : tension will decrease.

### 2. Adjustment of joint friction (fig.4)

Friction value is factory set, but can be adapted to the patient gait by screwing or unscrewing screw (F) with a 3mm hex wrench :

- clockwise : friction will increase
- anti-clockwise : friction will decrease.

A change can already be noticed with a 15° to 20° rotation of the hex wrench.

**Caution !** If the adjustment is too loose, it can lead to noises and a premature wear. On the contrary, a too tight adjustment can lead to unexpected braking, even without weight bearing.

### 3. Adjustment of brake sensitivity according to load (fig.5)

Brake sensitivity can be adjusted by screwing or unscrewing screw (S) with a 6 mm hex wrench :

- clockwise : a higher load is necessary to make the joint brake
- anti-clockwise : a lower load will make the joint brake.

Fig. 3


Fig. 4




Fig. 5

## **7 - Installation of the protective cover for the cosmetic (fig.6)**

Pull protective cover (**g**) on tube (**h**) and let it slide upwards until stud (**c1**) of mobile support (**c**) comes into hole (**g1**) of protective cover (**g**).

**Fig. 6**

## **8- Product life**

This component has been tested in accordance with standard ISO10328, load level P5 (100 kg, 220 lbs).

The cyclical test of 3 million cycles corresponds to a use of 4 or 5 years depending on the patient activity.

## **9- Maintenance**

A yearly control of the knee is recommended to check proper functioning and to clean it (dust). If necessary, proceed to a new adjustment of the knee.

If needed, extension assist 1M11510 can be ordered separately.

## **10 – Advice for use, maintenance and safety (inform the patient)**

In order not to damage the knee, do not use talcum powder, grease or a silicon spray to eliminate possible rubbing noise. Any use of these products may damage the mechanical parts, which may lead to function defects and risk of falling for the patient.

***The use of talcum powder, grease or silicon spray will render all claims against PROTEOR null and void.***

- It is strictly forbidden to tighten or loose any screw of this knee except for the 4 adjustment screws **R**, **F** and **S** (Fig.3, 4 and 5) that can be used by the prosthetist only.
- The warranty is not applicable to damages due to improper use, inadequate alignment, use in a very dusty environment and without adapted protection, or any inappropriate use.
- There is a risk of having a finger caught or clothing pinched in the knee. To avoid any risk of injury caused by the moving joint, take care that no one puts his finger near or inside the mechanism.
- It should be avoided to expose the knee in environments that can cause corrosion of metal parts (fresh water, sea water, chlorinated water, acids, etc.)
- In case of fall, shock or abnormal behaviour of the prosthesis, the knee should be checked by the prosthetist.
- The temperature range for knee use is from -10°C to +40°C (from 14°F to 104°F).

***Please give a copy of these instructions for use to the patient.***

## **11 – Recycling—End of life**

This product includes items of various materials : elastomer, titanium, steel and brass.

They must be recycled in compliance with the laws in force.

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