## CONECTOR FEMORAL ORIENTABLE 1K176

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### ASSEMBLING INSTRUCTIONS

#### I. CONTENTS


- 1 lamination anchor A with 4 hollow hexagonal grub screws M6x12 (A1)
- 1 turning base B with 4 hollow hexagonal grub screws M8x16 (B1)
- 4 hollow hexagonal grub screws M6x16 (A2).

#### II. SHAPING THE LAMINATION ANCHOR

- The prongs can be adapted to the form of the positive for aesthetic matter.
- CAUTION : Nevertheless when doing so take care not to distort the central hole because this could hinder rotation.

#### III. RECOMMENDATIONS FOR LAMINATION

- Remove screws A1 (M6x12) and replace them by screws A2 (M6x16). Let them protrude from about 4 mm. These screws will be used to find the position of the threaded holes when eliminating lamination excess.
- Stop the central hole with plastiline to avoid resin penetration.
- Laminate with your usual method (glass or carbon rovings through anchor holes and several layers of glass or carbon fabric inside and outside according to the resistance you wish).
- When resin setting is completed, remove lamination excess taking care not to damage the lower side and the threaded holes.
- Reset the 4 screws A1 (M6x12).


#### IV. CONNECTIONS

- This product can be installed on any "standard male pyramid" (e.g. 1K160).

#### V. ADJUSTMENTS

- Rotation of anchor A around turning base B can be locked by tightening the 4 screws A1 (advised tightening torque : 5 N.m).
- Act on the 4 screws B1 (advised tightening torque : 15 N.m) to adjust tilt.
- When adjustments are completed, use glue (ref. Proteor XC047) to fix them.

#### VI. USE

- Active patients involved in sports such as cycling or quick walking.
- Patients up to 125 kg.