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**ASSEMBLY INSTRUCTIONS  
THORACO-BRACHIAL BRACE KIT 6D050**


**I - USE**

The parts and accessories of this kit are used to connect the different sections of a thoraco-brachial brace previously shaped to the patient's dimensions on a positive mould.

The kit may be used for either right arm or left arm.


**II - KIT CONTENTS**

The kit includes four sub-assemblies which may also be ordered as replacement parts :  
6D050 = 6D05010 + 6D05020 + 6D05050 + 6D05040

### III - FASTENING THE KIT COMPONENTS TO THE «CORSET / ARM ORTHOSIS» ASSEMBLY

**Figure 1**


A - Temporarily fasten the corset, arm orthosis and forearm orthoses to the positive plaster mould.

B - Line up the shoulder brackets on the mould and shape them so that the pivot point goes through the **shoulder joint** axis marked on the mould (see figure 1).

Make sure the drilled parts of the shoulder brackets correctly match the shape of the corset and of the arm orthosis at the fastening points.

**Figure 2**


Fit the four fastening plates on the brackets in the intermediate setting position, then mark the position of the plates on the corset and arm orthosis with a felt-tip pen.

C - Line up the arm support on the arm and forearm orthoses so that the arm support axis matches the elbow **joint axis** marked on the positive mould (see figure 2).

Mark the position of the support on the arm and forearm orthoses with a felt-tip pen.

**Figure 3**


D - Remove the «corset / arm orthosis» assembly from the positive mould. Trim around the useful parts (taking care not to remove the marks of the previously positioned parts) and then fix in place **temporarily** :

a) the four shoulder bracket fastening plates with the eight M5x25 machine screws (see figure 3).


**Temporarily** refit the shoulder brackets with the eight HM5 nuts.

b) the arm support with the four 3x15 countersunk head screws and the four HM3 nuts (see figure 4).

**Figure 4**

c) the height adjustment strut to the arm support in one of the two holes (see figure 4) and to the corset with the two M5x25 machine screws and the two HM5 nuts (see figure 5).

**Figure 5**

#### IV - ADJUSTMENT


Test the «corset / arm orthosis» assembly on the patient. Modify the position of the arm support, or the attachment points or length of the height adjustment strut if necessary.

#### A - LENGTH OF THE HEIGHT ADJUSTMENT STRUT (see fig. 6 and 6 bis)


ALWAYS SUPPORT THE PATIENT'S ARM IN THE ORTHOSIS WHEN ADJUSTING THE STRUT .

1. Hold the arm and orthoses in the maximum raised position.
2. Unscrew the lower pivot (black knob) of the height adjustment strut
3. Extract the ringed rod from the sleeve and remove the black O-ring that acts as a position marker.
4. Replace the ringed rod on the hip support with its pivot (black knob).
5. Hold the ringed rod parallel to the sleeve and mark the groove that is opposite mark 12 (see figure 6).

**Figure 6**

**Figure 6bis**

6. Cut off the extra length at the next groove so as to keep sufficient insert length.
7. Put the black O-ring (see figure 6 bis) back in the groove marked previously.
8. Fit the ringed rod back in the sleeve and check that it slides properly.

## **B - STRAPS POSITIONNING**

Mark the positions of the straps on the orthoses :

- 25 mm strap on the forearm orthosis,
- 40 mm strap on the arm orthosis,
- pectoral strap on the corset.

## **V - FINISHING**

After making all the necessary settings and adjustments, replace the temporary fastenings as follows :

1. The four fastening plates : (x8) 3x16 copper rivets
2. The arm support (after shaping) : (x7) 3x16 copper rivets  
(Note : Sandwich the arm support between the arm orthosis and the height adjustment strut)
3. The support of the strut base : (x2) 4x20 copper rivets
4. The arm and forearm straps : (x4)  $\varnothing 4$  plain nickel rivets
5. The pectoral strap : (x2)  $\varnothing 4$  plain nickel rivets  
Caution : - Permanently tighten the screws and pivots of the shoulder joints  
- Permanently tighten the arm and support at the angle set by the practitioner  
- Fit the pads to the straps
6. Fasten the ventral strap (elastic backing with double tightening) on the corset.  
(Note : This strap does not come with the kit)

## **VI - PRECAUTIONS**

Check regularly that the pivots and retaining screws of the shoulder bracket are tight.

 ALWAYS SUPPORT THE PATIENT'S ARM IN THE ORTHOSIS WHEN ADJUSTING THE STRUT.

 AFTER ADJUSTING THE STRUT, MAKE SURE THE KNURLED SCREW IS LOCKED IN A GROOVE ON THE RINGED ROD, AND SECURE IT WITH PROTEOR XC046 ADHESIVE.

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