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**I PACKING CONTENT (Fig.1)**

- A lamination adapter (A) with its clamping screw M5, length 30 (A1) and its spheric washer (A2)
- A screw M5, length 40 (B) used with the dummy during lamination
- A dummy made of plastics (C)
- This instructions sheet..

II USE

This component is designed exclusively for orthopedic fitting of transfemoral amputees whose weight, carried load included, does not exceed 150 Kg.

III ADAPTATION TO THE SOCKET FORM

- The prongs of adapter 1K179 can be shaped with your usual method to match the socket form.
- Caution : While shaping the prongs of the connector, take care not to bend them with too low curvature radius and do not use forming tongs whose claws have sharp edges. To avoid distortion of the connector thread, we also recommend that adapter 1K179 be placed during shaping on a component with a male thread M36x1,5.

IV SOCKET LAMINATION (Fig.2)

- Remove screw A1 (M5x30) and replace it by screw B (M5x40)
- Fully tighten device C in the M36x1.5 tapped part
- Protect slot F with plastiline
- Before placing the adapter on the plaster model, cover the latter with several layers of glass or carbon fabrics
- Position the adapter and check its alignment very carefully, then put some layers of glass or carbon fabrics over it. Tie off these fabrics layers at the thread base
- Perform lamination with your usual method
- After resin polymerization, remove the excess material to release device C and screw B
- Remove screw B and device C
- Remove any trace of plastiline.
- Place screw B, tighten it fully, then unscrew it completely. Repeat this several times to thoroughly remove any resin residue from the tapped part.
- Then install screw A1 (M5x30).

V CONNECTION (Fig.3)

Adapter 1K179 can be connected to any component D with a male thread M36x1,5 (1K183, 1K155, ...).

Caution : never unscrew the lamination adapter more than a quarter of a turn and check that it is screwed by 7 thread steps at least. Otherwise the component would get weak.

Tighten screw A1 up to a 7.4 Nm torque.

VI CHECK

Functions of adapter 1K179 are to be checked in regular intervals and at least once a year.

VII RECYCLING

The main part of this component is made of stainless steel and the screws are made of steel. Each of these items should be recycled according to the law in force.

Fig.1


Fig.2



Fig.3

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