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

**1N106 AND 1N107 CONNECTORS FOR
1N06 AND 1N07 THERMOFORMED SOCKETS :
Instructions for Use**

A) CONTENTS (Fig.1)

The system includes :

- x1 spherical base **1** with its screws **1a and 1b**
- x1 nut washer **2** with its screws **2a**
- x1 milled head screw **3**, 25 mm length
- x1 milled head screw **4**, 30 mm length
- x2 headless hexagonal socket screws **5** (see Use §E => Fig. 5b) :
M8 x 20 for 1N106
M6 x 16 for 1N107

B) USE

These connectors are used for patients whose maximum weight, carried load included, does not exceed :


- 100 kg for **1N106**
- 45 kg for **1N107**.

C) CONNECTIONS (Fig.2)

These connectors are designed to :

- receive between spherical base **1** and nut washer **2 a** thermoformed socket **A** of the PROTEOR range :
- **1N06**** for connector **1N106**
- **1N07**** for connector **1N107**.
- be connected through spherical base **1** to any standard male pyramid connector **B** :
- of the adult range for 1N106 (1K163, ect...)
- of the child range for 1N107 (1K164, ect...)

**** CAUTION :** Using a thermoformed socket in plastics outside the PROTEOR range, or any other type of socket (even from the PROTEOR range), may weaken the prosthesis and induce a risk for the patient.

Depending on the socket thickness, select screws **3** or **4** supplied to assemble the prosthesis :


- For a socket bottom whose thickness (X) is up to 9 mm, select screw **3** supplied (25 mm length)
- For a socket bottom whose thickness (X) ranges from 9 to 15 mm, select screw **4** supplied (30 mm length).

D) SOCKET TRANSLATION ADJUSTMENT (Fig.3 / 4)

1N106 and 1N107 connectors allow to adjust the socket translation within a range of ± 5 mm.

Before starting adjustments, screws 1a of base 1 and screws 2a of washer nut 2 must be so undone that their tips do not get in touch with socket A.


To adjust socket translation, central screw 3 or 4 must be slightly undone to allow socket A to slide between base 1 and screw nut 2.

Once the adjustment is completed, slightly screw central screw 3 or 4 to lock the assembly, then tighten screws 1a and 2a using a 2 mm hex wrench, so that their tips get into the plastics of socket A (see fig. 4).

To make this adjustment reliable, tighten central screw 3 or 4 :

- until you reach a 22 Nm torque, using a 6 mm hex wrench for 1N106
- until you reach a 17 Nm torque, using a 5 mm hex wrench for 1N107

E) SOCKET ANGULATION ADJUSTMENT (Fig.5a / 5b)


To adjust angulation of socket A, turn screws 1b of base 1 using an hex wrench (4 mm for 1N106 or 3 mm for 1N107) until you reach the desired adjustment.

If one of the 1b screws is too short, replace it by one of the longer screws 5 included in the pack (M8x20 for 1N106 and M6x16 for 1N107)

See Fig. 5b

Apply adhesive XC047 on the setting screws of the pyramid, then tighten each screw up to the following torque :

- 15 Nm for 1N106
- 7,5 Nm for 1N107.

F) RECYCLING

- Spherical base 1 and washer nut 2 are made of aluminum alloy.
- Screws 1a, 1b, 2a, 3, 4 and 5 are made of zinc-coated steel.

Each of them should be recycled according to the law in force.

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