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

1- Description and functions

The **Dynatrek** foot includes the following main components :

- An upper carbon fiber blade that allows energy storage and return, as well as inversion/eversion at forefoot.
- A carbon fiber plantar keel that ensures a progressive and smooth step rollover.
- Two elastomer wedges for heel stiffness adjustment.
- A polyurethane foot shell (1A60020) which transmits the mechanical stresses between the shoe and the composite structure.
- A technical fiber sock limiting noises and wear and tear between the foot shell and the structure.
- A malleolar clip that ensures a good assembly with the cosmetic foam.

Thanks to the combination of these components, the foot has a comfortable and flexible behaviour.

2- Indications, maximum weight and connecting components

The **Dynatrek** foot is designed exclusively for lower limb prosthesis.

It is appropriate for patients with normal or high activity level whose weight (carried load included) complies with the table opposite.

Module selection guide						
	Patient weight (Kg)					
Activity	45 to 59	60 to 74	75 to 89	90 to 104	105 to 124	125 to 150
Normal	1	2	3	4	5	6
High	2	3	4	5	6	-

The authorized maximum load is 150Kg (330 lbs). This module has been tested according to standard NF EN ISO 10328-2006 for load level P7.

The foot is resistant to fresh, sea or chlorinated water provided that the recommendations of § 9 are respected.

The **Dynatrek** foot is designed for a 10 mm heel height.

It includes a male pyramid that allows the use of pyramid receivers (1D41 or 1D41-HD). Thanks to these connectors, the heel height can be modified from 0 to 20 mm.

3- Assembling (without the patient)

Place under the foot a wedge of the same height as the shoe heel. Fix the pyramid receiver onto the foot, taking care that its axis be perpendicular to the ground (fig. 1).


Fig. 1

4- Alignment (with the patient)

Place the foot so that the load line goes through the posterior third of the foot (fig. 2).



Fig. 2

- Build height :**
- 154mm from T22 to T25
 - 164mm from T26 to T30

- Increasing equinus provides a more dynamic alignment. It increases energy return at step completion.
- If the foot seems too long during transfer to the forefoot or if heel support is too flexible, you have to position the foot behind the socket.
- If heel support seems to hard, or if the foot seems too short during transfer to the forefoot, you have to position the foot ahead of the socket.

5- Heel adjustment

The foot response at heel strike can be modified using the included elastomer wedges. The hardness of each wedge will provide more or less firm support depending on the patient's need : red = moderate hardness, yellow = high hardness. Carry out tests by scotching a wedge in the gap between the two blades. Once the correct hardness is determined, fix definitively the wedge with cyanoacrylate glue (Ref. Proteor XC050 / Sicomet 8300) following the steps below :

6- Finishing : To avoid rubbing noises, pull the sock up to the basis of the pyramid and fix it with a fine string for example.

The malleolar clips shown in the following table ensure connection to the cosmetic shin foam and allow to remove the foot easily and rapidly.

	Foot size (cm)	22 to 25	26 to 27	28 to 30
	Product Nr	1A40013	1A40014	1A40015

Installation : Cut the cosmetic foam to the desired length, position the foam on the foot and draw the malleolus contour. Inside the foam mill a space that will receive the selected connecting part and the foot part protruding from the shell, and at the same time take care that thickness be sufficient to allow for sticking. Position the malleolar clip carefully and clip it into the cosmetic foam. Clean contact surfaces with isopropyl alcohol and stick the malleolar clip on the foam with XC090 neoprene glue. After drying, finish the outer shape as usually.

7- Product life :

This component has been tested in accordance with standard NF EN ISO 10328, the cyclical test of 3 million footsteps corresponds to a use of 2 to 3 years depending on the patient activity.

An annual check of the knee by a CPO is recommended.


8- Maintenance :

The foot shell, the sock, the malleolar clip and the heel wedge are the only foot components that may need to be replaced (if worn out, lost, ...)

Foot shell :


The foot shell can be removed and installed again using 1A40019 shoe-horn.

Use exclusively the spare kit 1A60020-D22 including the shell, the sock and the malleolar clip.

Wedge :

Use wedge 1A60014 (red = hardness 40ShA) or 1A60015 (yellow = hardness 60ShA) depending on the level of comfort required at the heel.

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9– Advice for use, maintenance and safety

Clean the foot shell with a slightly damp cloth or sponge.

Use temperature : - 20°C / + 60°C.

No maintenance operation such as lubrication, work on the screws or other parts is required for this foot except for the operations indicated in §8 that must be performed by the CPO.

A use of the foot that does not adhere to your CPO's recommendations can damage foot components (carrying heavy loads for example).

In case of an abnormal behaviour of your foot, or if you feel changes in your foot features, or if the foot sustained a hard shock, contact your CPO.

The presence of sand, particles or any foreign body inside the foot could alter its behaviour and damage it.

After use in water or in a very dirty environment :

Remove the foot shell and take off the sock.

Rinse the foot and the other parts thoroughly with clear water to eliminate particles (mud, sand, ...).

Dry the different items before replacing the sock and the shell.

10– Recycling

The different items of the foot are special wastes and must be recycled according to the laws in force.

11– Regulatory information for France

Testing laboratory Reports—Dates	Nr of points	Permanent deformation (mm)	Inversion / Eversion angle	Sagittal flexion amplitude
C.E.R.A.H Test reports : N° 16-005-A of May 2nd, 2016	133,32	Heel : 1,17mm Forefoot : 1,43mm	Inversion : 5,04° Eversion : 8,53°	Plantar flexion : 7,98° Dorsiflexion : 8,52°

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