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

**1 - Description and functions**

The **DynaCity** foot includes the following main components :

- An upper carbon fiber blade that allows energy return.
- A carbon fiber plantar blade that ensures a progressive and smooth step rollover.
- An intermediate carbon fiber blade allowing to get flat foot very well.
- The 3 blades are split to provide good inversion-eversion at forefoot.
- A polyurethane foot shell (1A60020) which transmits the mechanical stresses between the shoe and the composite structure.
- A technical fiber sock limiting noises and wear and tear between the foot shell and the structure.
- A malleolar clip that ensures a good assembly with the cosmetic foam.

Thanks to the combination of these components, the foot has a comfortable and flexible behaviour.

**2- Indications, maximum weight and connecting components**

The **DynaCity** foot is designed exclusively for lower limb prosthesis.

It is appropriate for patients with low or moderate activity level whose weight (carried load included) complies with the table opposite. The maximum load has been tested according to standard ISO-22675.

<b>Category Selection Guide</b>					
<b>Size</b>	<b>Patient weight (Kg)</b>				
	45 to 59	60 to 74	75 to 90	91 to 107	108 to 125
<b>22</b>	1	1	2	2	NA
<b>23</b>	1	1	2	2	NA
<b>24</b>	1	1	2	2	NA
<b>25</b>	2	2	3	3	4
<b>26</b>	2	2	3	3	4
<b>27</b>	2	3	3	4	4
<b>28</b>	2	3	3	4	4
<b>29</b>	4	4	4	5	5
<b>30</b>	4	4	4	5	5

The **DynaCity** foot is designed for a 10 mm heel height.


The foot is resistant to fresh water only, provided that the recommendations of § 8 are respected.

\* Categories are not modules : not only the patient weight, but also the foot size must be taken into account to determine the category.

It includes a male pyramid that allows the use of standard adapters. Thanks to these adapters, the heel height can be modified from 0 to 20 mm.



**3- Assembling (without the patient)**

Place under the foot a wedge of the same height as the shoe heel. Fix the pyramid receiver onto the foot, taking care that its axis be perpendicular to the ground (fig. 1).

**4- Alignment (with the patient)**

Place the foot so that the load line goes through the posterior third of the foot (fig. 2).

- Increasing equinus provides a more dynamic alignment. It increases energy return at step completion.
- If the foot seems too long during transfer to the forefoot or if heel support is too flexible, you have to position the foot behind the socket.
- If heel support seems too hard, or if the foot seems too short during transfer to the forefoot, you have to position the foot ahead of the socket.

**5 – Finishing :** To avoid rubbing noises, pull the sock up to the basis of the pyramid and fix it with a fine string for example.

The malleolar clips shown in the following table ensure connection to the cosmetic shin foam and allow to remove the foot easily and rapidly.

	Foot size (cm)	22 to 25	26 to 27	28 to 30
	Product Nr	1A40013	1A40014	1A40015

**Installation :** Cut the cosmetic foam to the desired length, position the foam on the foot and draw the malleolus contour. Inside the foam mill a space that will receive the selected connecting part and the foot part protruding from the shell, and at the same time take care that thickness be sufficient to allow for sticking. Position the malleolar clip carefully and clip it into the cosmetic foam. Clean contact surfaces with isopropyl alcohol and stick the malleolar clip on the foam with XC090 neoprene glue. After drying, finish the outer shape as usually.

### **6– Product life :**

This component has been tested in accordance with standard ISO 22675, the cyclical test of 2 million cycles corresponds to a use of 2 to 3 years depending on the patient activity.

An annual check of the knee by a CPO is recommended.



### **7- Maintenance :**

The foot shell, the sock and the malleolar clip are the only foot components that may need to be replaced (if worn out, lost, ...)

#### Foot shell :

The foot shell can be removed and installed again using 1A40019 shoe-horn.

Use exclusively the spare kit 1A60020-D22 including the shell, the sock and the malleolar clip.

### **8– Advice for use, maintenance and safety**

Clean the foot shell with a slightly damp cloth or sponge.

Use temperature : - 20°C / + 60°C.

No maintenance operation such as lubrication, work on the screws or other parts is required for this foot except for the operations indicated in §7 that must be performed by the CPO.

A use of the foot that does not adhere to your CPO's recommendations can damage foot components (carrying heavy loads for example).

In case of an abnormal behaviour of your foot, or if you feel changes in your foot features, or if the foot sustained a hard shock, contact your CPO.

The presence of sand, particles or any foreign body inside the foot could alter its behaviour and damage it.

After use in water or in a very dusty environment :

Remove the foot shell and take off the sock.

Rinse the foot and the other parts thoroughly with clear water to eliminate particles (mud, sand, ...).

Dry the different items before replacing the sock and the shell.

## **9– Recycling**

The different items of the foot are special wastes and must be recycled according to the laws in force.

## **10– Regulatory information for France**

Testing laboratory Reports—Dates	Nr of points	Permanent deformation (mm)	Inversion / Eversion angle	Sagittal flexion amplitude
C.E.R.A.H Test reports : N° 16-162-B of March 14th, 2017	117,89	Heel : 1,36 mm Forefoot : 0,81mm	Inversion : 6,39° Eversion : 5,92°	Plantar flexion : 8,11° Dorsiflexion : 10,26°

**PROTEOR**  
6 rue de la Redoute - 21850 Saint-Apollinaire - France  
Tél. 33 (0) 3 80 78 42 10 - Fax 33 (0) 3 80 78 42 15  
ht.orthopedie@proteor.com - www.proteor.com

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