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

Abb. III


## INSTRUCTIONS FOR FOUR-AXIS KNEE JOINT WITH STABILITY ADJUSTMENT 1M03

### JOINT FRICTION ADJUSTMENT (Fig I)

Screw in or out with a hex wrench to achieve the desired friction. Please apply the same amount of adjustment to each screw.

Fig. I

### EXTENSION BEARING ADJUSTMENT (Fig II)

The knee joint is delivered with two extension bearings : the regular one is assembled with knee joint, the hard one could be used instead of the regular one.

Screw in or out the bearing screw with a hex wrench ; Make that adjustment when the patient weight is applied to the prosthesis in order to keep the bearing correctly in its housing.

Fig. II


- 1.Unlock the central screw with a 2 mm hex wrench.
2. Tighten or loosen the main screw with a 4 mm hex wrench.
3. Lock the central screw (2.5 Nm) with a 2 mm hex wrench.

### STABILITY ADJUSTMENT (Fig II)

Stability of knee 1M03 is adjusted by moving forward or backward axis 4 by means of piece 3 which is fastened to tube-holder 5 with screw 1.

The highest stability is achieved when heads of screws 1 and 2 are located on the same side of the knee-joint (as shown opposite).

The regular stability is achieved when heads of screws 1 and 2 are located on opposite sides of the knee joint.

### EXTENSION SPRING ADJUSTMENT (Fig II)

Tighten or loosen the screw with a flat blade screwdriver to get the desired extension force.

**Caution** ! Then check that the knee can achieve maximal flexion. Otherwise unscrew slightly the screw.

### ASSEMBLING—ALIGNMENT (Fig III)


The position of the knee axis relative to the prosthesis load line is paramount for the joint to work properly and to last a long time. We recommend the alignment shown opposite.

### USE

Slow-walking patients. Sedentary patients. Low activity level.

Patients (carried load included) up to 100 Kg.

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