**1M01 - 1M10
SINGLE AXIS LOCKING KNEE
INSTRUCTIONS FOR USE**

**1M01 - 1M10
GENOU MONOAXIAL A VERROU
NOTICE DE MONTAGE**

USE

Slow-walking or sedentary patients with low activity level.
Patients up to 100 kg (carried load included).
Patients requiring the maximum safety provided by a locked knee.

ALIGNMENT (Fig. 2)

The alignment shown opposite is highly recommended.

BEARING ADJUSTMENT (Fig. 1)

Tighten or loosen screw (1) with a 3 mm hex wrench until you achieve the desired adjustment.

CAUTION ! Check systematically that the knee locks in extension (locking noise and free motion range of the lever = about 2 mm)

CABLE MOUNTING (Fig. 3)

- Place cross piece (3) in the holes of lever (2).
- Insert screw (4) in cross piece (3).
- Fit nut (5) on screw (4).
- Insert rivet (6) in the hole of nut (5).
- Put cable (7) in rivet (6).
- Adjust the length of cable (7).
- Tighten screw (4) to clamp the cable
- Cut off the excess material.
- Place cap (8) on nut (5).

CAUTION : In maximal flexion the socket must necessarily get in contact with the collar of the tube-holder. Any other contact is forbidden.

FRICITION ADJUSTMENT (Fig. 1)

1. Unlock central screw (9) with a 2 mm hex wrench.
2. Tighten or loosen main screw (10) with a 4 mm hex wrench.
3. Lock central screw (9) : tightening torque 2.5 Nm.

HOW TO CONVERT 1M01 → 1M10 (Fig. 4)

Remove two screws (11) and take care not to lose the joint axes.
Select position A or B for handle (12).
Fit covers (13) on handle (12).
Insert the ends of handle (12) in the holes of lever (2).
Connect them together with screws (14).

HOW TO CONVERT 1M1 → 1M01 (Fig. 4)

Remove screws (14), covers (13) and handle (12).
Mount the cable as indicated above.

Replace screws (14) by hex headless screws provided with the cable !

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