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

- Active patients.
- Patients up to 125 kg.

II. CONTENTS

- A lamination connector **A** with 4 hollow hexagonal headless screws M6 x12 (**A1**)
- A turning male pyramid **B**
- 4 hollow hexagonal headless screws M6x16 (**A2**)

III. SHAPING THE LAMINATION ANCHOR

- The prongs can be adapted to the form of the positive to improve the reliability of the link between the connector and the socket and for aesthetic matter. Caution : Nevertheless when doing so take care not to distort the central hole, because this could hinder rotation.

III. RECOMMENDATIONS FOR LAMINATION

- Remove screws **A1** (M6x12) and replace them by screws **A2** (M6x16). Let them protrude from about 4 mm. These screws will be used to find the position of the threaded holes when eliminating resin excess.
- Stop the central hole with plastiline to avoid resin penetration.
- Laminate as usual (glass or carbon rovings through connector holes and several layers of glass or carbon fabric inside and outside according to the resistance you wish).
- After polymerization and cooling, remove resin excess taking care not to damage lower side **C** and threaded holes.
- Reset the 4 screws **A1** (M6x12)

IV. ADJUSTMENTS

- Rotation of connector **A** around turning base **B** can be locked by tightening the 4 screws **A1** (advised tightening torque : 5Nm)
- We recommend to secure the finalized adjustments with glue (Proteor ref. XC047).

V. CONNECTIONS

- This product is compatible with any "standard female pyramid".